



05.0017.10

GB FOOTSWITCH FOR SUCTION PUMP
Vacuum gauge regulation.

F PEDALE POUR ASPIRATEUR
Régulation du vide.

SUITABLE FOR / UTILISABLE AVEC

F-30.XX/P
F-31.XX/P
F-35.XX/P
F-X50 Series / Série
F-60 Series / Série
F-90 Series / Série
F-100 Series / Série



JARS

GB JARS
Autoclavable jars for suction pumps made of polycarbonate (Makrolon) autoclavable at 121°C. Silk-screened graduated scale resistant to aggressive liquids and disinfectants.

F BOCAUX
Bocaux pour aspirateurs. Polycarbonate (Makrolon) autoclavable à 121°C. Graduation résistante aux liquides agressifs et désinfectants.

TECHNICAL SPECIFICATIONS / SPECIFICATIONS TECHNIQUES

AUTOCLAVABLE 121°	SINGLE USE/USAGE UNIQUE
05.0001 1 L	D1LT 1 L
05.0002 2 L	D2LT 2 L
05.0004 4 L	D4LT 4 L



M043002 (1L) / M045002 (2L)

GB MEDIFLOW BAG:
Non-sterile bag for collecting body fluids in surgery disposable. Medical latex free polyethylene. Minimum pressure -560 mmHg (-750 mbar). Avoid contact with the aspirated secretions.

F MEDIFLOW BAG
Sac non stérile de collecte des liquides corporels en chirurgie, à usage unique. Polyéthylène médical sans latex. Pression minimale -560 mmHg (-750 mbar). Evite tout contact avec les sécrétions aspirées.



SILICONE TUBES

GB SILICONE TUBES
Silicone tubes for all suction pumps

F TUBES EN SILICONE
Tubes en silicone pour tous les aspirateurs.

REFERENCES

SP.27.01	Silicone tube / Tube en silicone Ø 11 x 6mm, 10 mt. Shore 75
SP.27.02	Silicone tube / Tube en silicone Ø 14 x 8mm, 10 mt. Shore 60
SP.27.03	Silicone tube / Tube en silicone Ø 13 x 7mm, 10 mt. Shore 60

Available different sizes on demande
Dimensions différentes disponibles sur demande



05.0016

GB ANTIBACTERIAL FILTER
Disposable.

F FILTRE ANTIBACTÉRIEN
Jetable.